

# Anti-Pentraxin 3/PTX3 [R05-3Q6] Monoclonal Antibody

## AGMB02660

### Description

---

This Anti-Pentraxin 3/PTX3 [R05-3Q6] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB02660
<b>Contents:</b>	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	PTX3, TNFAIP5, TSG14, Pentraxin-related protein PTX3, Pentaxin-related protein PTX3, Tumor necrosis factor alpha-induced protein 5, TNF alpha-induced protein 5, Tumor necrosis factor-inducible gene 14 protein, TSG-14
<b>Applications:</b>	<b>WB</b> <b>FC</b>
<b>Research Area:</b>	Cardiovascular
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Human, Mouse
<b>Clone:</b>	R05-3Q6
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P26022
<b>Immunogen:</b>	A synthesized peptide derived from human Pentraxin 3
<b>Gene ID:</b>	5806

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** PTX3  
**Host Species:** Rabbit  
**Isotype:** IgG  
**Purification:** Affinity Chromatography  
**Conjugated:** Unconjugated  
**Modification:** Unmodified  
**Molecular Weight:** Calculated MW: 42 kDa, Observed MW: 42 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

	Application	Antibody Dilution Ratio
<b>Antibody Dilution Ratio:</b>	WB	1:1000-1:2000
	FC	1:50-1:100

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.