

Anti-PDE4D [R08-9G2] Monoclonal Antibody

AGMB02915

Description

This Anti-PDE4D [R08-9G2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB02915
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	cAMP-specific 3', 5' -cyclic phosphodiesterase 4D, DPDE3, PDE43, PDE4D, DPDE3, PDE4A, PDE4B, PDE4C,
Applications:	
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R08-9G2
Clonality:	Monoclonal Antibody
SwissProt ID:	Q08499
Immunogen:	A synthesized peptide derived from human PDE4D
Gene ID:	5144

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: PDE4D

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 91 kDa, Observed MW: 91 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.