

Anti-Ctip1/BCL-11A [R07-2N7] Monoclonal Antibody

AGMB02936

Description

This Anti-Ctip1/BCL-11A [R07-2N7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB02936
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	EVI9, CTIP1, ZNF856, HBFQTL5, BCL11A-L, BCL11A-S, BCL11a-M, BCL11A-XL, 2810047E18Rik, B cell CLL/lymphoma 11A (zinc finger protein), B cell CLL/lymphoma 11A (zinc finger protein) isoform 2,
Applications:	WB IHC-P ICC/IF FC
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R07-2N7
Clonality:	Monoclonal Antibody
SwissProt ID:	Q9H165
Immunogen:	A synthesized peptide derived from human Ctip1
Gene ID:	53335

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: BCL11A

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 91 kDa, Observed MW: 120 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IHC	1:100-1:200
	IF	1:50-1:200
	FC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.