

# Anti-LMAN1 [R01-8R1] Monoclonal Antibody

## AGMB02947

### Description

---

This Anti-LMAN1 [R01-8R1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |                                                                                                                                                                              |
|-----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>SKU:</b>           | AGMB02947                                                                                                                                                                    |
| <b>Contents:</b>      | 20µl, 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml)                                                                                                                          |
| <b>Synonyms:</b>      | ERGIC5, F5F8D, FMFD1, Gp58, Lman1, MCFD1, MR60, Endoplasmic reticulum golgi intermediate compartment protein 53, ER-Golgi intermediate compartment 53 kDa protein, ERGIC-53, |
| <b>Applications:</b>  | <b>WB</b> <b>IHC-P</b> <b>ICC/IF</b>                                                                                                                                         |
| <b>Research Area:</b> | Signal Transduction                                                                                                                                                          |
| <b>Form:</b>          | Liquid                                                                                                                                                                       |

### Antibody Data

---

|                      |                                                |
|----------------------|------------------------------------------------|
| <b>Reactivity:</b>   | Human, Mouse, Rat                              |
| <b>Clone:</b>        | R01-8R1                                        |
| <b>Clonality:</b>    | Monoclonal Antibody                            |
| <b>SwissProt ID:</b> | P49257                                         |
| <b>Immunogen:</b>    | A synthesized peptide derived from human LMAN1 |
| <b>Gene ID:</b>      | 3998                                           |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

|                          |                                            |
|--------------------------|--------------------------------------------|
| <b>Gene Name:</b>        | LMAN1                                      |
| <b>Host Species:</b>     | Rabbit                                     |
| <b>Isotype:</b>          | IgG                                        |
| <b>Purification:</b>     | Affinity Chromatography                    |
| <b>Conjugated:</b>       | Unconjugated                               |
| <b>Modification:</b>     | Unmodified                                 |
| <b>Molecular Weight:</b> | Calculated MW: 58 kDa, Observed MW: 58 kDa |

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

| <b>Antibody Dilution Ratio:</b>           | Application                                                                                                                                                                                                                                             | Antibody Dilution Ratio |
|-------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------|
|                                           | WB                                                                                                                                                                                                                                                      | 1:1000-1:2000           |
|                                           | IHC                                                                                                                                                                                                                                                     | 1:100-1:200             |
|                                           | IF                                                                                                                                                                                                                                                      | 1:50-1:200              |
| <b>Protein Quantification (Optional):</b> | To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <a href="https://www.assaygenie.com/bradford-protein-assay-protocol/">https://www.assaygenie.com/bradford-protein-assay-protocol/</a> to view the full protocol. |                         |