

# Anti-Parathyroid Hormone Receptor 1 [R05-7Y8] Monoclonal Antibody

## AGMB02951

### Description

---

This Anti-Parathyroid Hormone Receptor 1 [R05-7Y8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB02951
<b>Contents:</b>	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	PTH1R, PTHR, PTHR1, Parathyroid hormone/parathyroid hormone-related peptide receptor, PTH/PTHrP type I receptor, PTH/PTHr receptor, Parathyroid hormone 1 receptor, PTH1 receptor
<b>Applications:</b>	<b>WB</b> <b>IP</b>
<b>Research Area:</b>	Signal Transduction
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Human
<b>Clone:</b>	R05-7Y8
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	Q03431
<b>Immunogen:</b>	A synthesized peptide derived from human Parathyroid Hormone Receptor 1

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene ID:** 5745

**Gene Name:** PTH1R

**Host Species:** Rabbit

**Isotype:** IgG

**Purification:** Affinity Chromatography

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** Calculated MW: 66 kDa, Observed MW: 66 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

**Antibody Dilution Ratio:**

Application	Antibody Dilution Ratio
WB	1:1000-1:2000
IP	1:50

**Protein Quantification (Optional):**

To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.