

# Anti-VPS26 [R09-6G2] Monoclonal Antibody

## AGMB02976

### Description

---

This Anti-VPS26 [R09-6G2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |                                                     |
|-----------------------|-----------------------------------------------------|
| <b>SKU:</b>           | AGMB02976                                           |
| <b>Contents:</b>      | 20µl, 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml) |
| <b>Synonyms:</b>      | FLJ12930, HB 58, HB58, VPS26                        |
| <b>Applications:</b>  | <b>WB</b> <b>ICC/IF</b> <b>IP</b>                   |
| <b>Research Area:</b> | Signal Transduction                                 |
| <b>Form:</b>          | Liquid                                              |

### Antibody Data

---

|                      |                                                |
|----------------------|------------------------------------------------|
| <b>Reactivity:</b>   | Human, Mouse, Rat                              |
| <b>Clone:</b>        | R09-6G2                                        |
| <b>Clonality:</b>    | Monoclonal Antibody                            |
| <b>SwissProt ID:</b> | O75436                                         |
| <b>Immunogen:</b>    | A synthesized peptide derived from human VPS26 |
| <b>Gene ID:</b>      | 9559                                           |
| <b>Gene Name:</b>    | VPS26A                                         |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

|                          |                                            |
|--------------------------|--------------------------------------------|
| <b>Host Species:</b>     | Rabbit                                     |
| <b>Isotype:</b>          | IgG                                        |
| <b>Purification:</b>     | Affinity Chromatography                    |
| <b>Conjugated:</b>       | Unconjugated                               |
| <b>Modification:</b>     | Unmodified                                 |
| <b>Molecular Weight:</b> | Calculated MW: 38 kDa, Observed MW: 38 kDa |

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

| <b>Antibody Dilution Ratio:</b> | Application | Antibody Dilution Ratio |
|---------------------------------|-------------|-------------------------|
|                                 | WB          | 1:1000-1:2000           |
|                                 | IF          | 1:50-1:200              |
|                                 | IP          | 1:50                    |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.