

# Anti-ILF3 [R02-5S3] Monoclonal Antibody

AGMB03058

## Description

---

This Anti-ILF3 [R02-5S3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

**SKU:** AGMB03058

**Contents:** 20µl, 50µl, 100µl  
Bradford Reagent: 1 vial (2ml)

**Synonyms:** ILF3, DRBF, MPHOSPH4, NF90, Interleukin enhancer-binding factor 3, Double-stranded RNA-binding protein 76, DRBP76, M-phase phosphoprotein 4, MPP4, Nuclear factor associated with dsRNA, NFAR, Nuclear factor of activated T-cells 90 kDa, NF-AT-90, Translational control protein 80, TCP80

**Applications:**  WB IHC-P ICC/IF FC IP

**Research Area:** Signal Transduction

**Form:** Liquid

## Antibody Data

---

**Reactivity:** Human, Mouse

**Clone:** R02-5S3

**Clonality:** Monoclonal Antibody

**SwissProt ID:** Q12906

**Immunogen:** A synthesized peptide derived from human ILF3

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene ID:** 3609

**Gene Name:** ILF3

**Host Species:** Rabbit

**Isotype:** IgG

**Purification:** Affinity Chromatography

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** Calculated MW: 95 kDa, Observed MW: 95 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

<b>Antibody Dilution Ratio:</b>	<b>Application</b>	<b>Antibody Dilution Ratio</b>
	WB	1:1000-1:2000
	IHC	1:100-1:200
	IF	1:50-1:200
	FC	1:50-1:100
	IP	1:50

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.