

# Anti-CLASP1 [R01-6Z6] Monoclonal Antibody

## AGMB03102

### Description

---

This Anti-CLASP1 [R01-6Z6] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |                                                                                                                                                           |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>SKU:</b>           | AGMB03102                                                                                                                                                 |
| <b>Contents:</b>      | 20µl, 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml)                                                                                                       |
| <b>Synonyms:</b>      | CLASP1, KIAA0622, MAST1, CLIP-associating protein 1, Cytoplasmic linker-associated protein 1, Multiple asters homolog 1, Protein Orbit homolog 1, hOrbit1 |
| <b>Applications:</b>  | <b>WB</b> <b>IHC-P</b> <b>ICC/IF</b> <b>FC</b>                                                                                                            |
| <b>Research Area:</b> | Cell Biology                                                                                                                                              |
| <b>Form:</b>          | Liquid                                                                                                                                                    |

### Antibody Data

---

|                      |                                                 |
|----------------------|-------------------------------------------------|
| <b>Reactivity:</b>   | Human, Mouse, Rat                               |
| <b>Clone:</b>        | R01-6Z6                                         |
| <b>Clonality:</b>    | Monoclonal Antibody                             |
| <b>SwissProt ID:</b> | Q7Z460                                          |
| <b>Immunogen:</b>    | A synthesized peptide derived from human CLASP1 |
| <b>Gene ID:</b>      | 23332                                           |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** CLASP1

**Host Species:** Rabbit

**Isotype:** IgG

**Purification:** Affinity Chromatography

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** Calculated MW: 169 kDa, Observed MW: 169 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

| <b>Antibody Dilution Ratio:</b> | <b>Application</b> | <b>Antibody Dilution Ratio</b> |
|---------------------------------|--------------------|--------------------------------|
|                                 | WB                 | 1:1000-1:2000                  |
|                                 | IHC                | 1:100-1:200                    |
|                                 | IF                 | 1:50-1:200                     |
|                                 | FC                 | 1:50-1:100                     |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.