

# Anti-RAIDD [R04-9P3] Monoclonal Antibody

## AGMB03114

### Description

---

This Anti-RAIDD [R04-9P3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB03114
<b>Contents:</b>	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	CRADD, MGC9163, RAIDD, Death adaptor molecule RAIDD, Death domain containing protein CRADD, CASP2 and RIPK1 domain containing adaptor with death domain, Caspase and RIP adapter with death domain, Caspase and RIP adaptor with death domain,
<b>Applications:</b>	<b>WB</b> <b>ICC/IF</b> <b>FC</b> <b>IP</b>
<b>Research Area:</b>	Cell Biology
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Human, Rat
<b>Clone:</b>	R04-9P3
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P78560
<b>Immunogen:</b>	A synthesized peptide derived from human RAIDD

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene ID:** 8738

**Gene Name:** CRADD

**Host Species:** Rabbit

**Isotype:** IgG

**Purification:** Affinity Chromatography

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** Calculated MW: 23 kDa, Observed MW: 23 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

<b>Antibody Dilution Ratio:</b>	<b>Application</b>	<b>Antibody Dilution Ratio</b>
	WB	1:1000-1:2000
	IF	1:50-1:200
	FC	1:50-1:100
	IP	1:50

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.