

Anti-Peroxiredoxin 5 [R09-1K6] Monoclonal Antibody

AGMB03166

Description

This Anti-Peroxiredoxin 5 [R09-1K6] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03166
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	ACR1, AOEB166, B166, Peroxiredoxin V, Peroxiredoxin-5, PLP, PMP20, PRDX5, PRDX6, Prx-V, PRXV, SBB110, Thioredoxin reductase
Applications:	WB IHC-P FC IP
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R09-1K6
Clonality:	Monoclonal Antibody
SwissProt ID:	P30044
Immunogen:	A synthesized peptide derived from human Peroxiredoxin 5
Gene ID:	25824

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: PRDX5

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 22 kDa, Observed MW: 17 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IHC	1:100-1:200
	FC	1:50-1:100
	IP	1:50

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.