

Anti-PDLIM1/CLP36 [R07-5G4] Monoclonal Antibody

AGMB03255

Description

This Anti-PDLIM1/CLP36 [R07-5G4] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

| | |
|-----------------------|---|
| SKU: | AGMB03255 |
| Contents: | 20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml) |
| Synonyms: | PDLIM1, CLIM1, CLP36, PDZ and LIM domain protein 1, C-terminal LIM domain protein 1, Elfin, LIM domain protein CLP-36 |
| Applications: | WB ICC/IF FC IP |
| Research Area: | Epigenetics |
| Form: | Liquid |

Antibody Data

| | |
|----------------------|---|
| Reactivity: | Human, Mouse, Rat |
| Clone: | R07-5G4 |
| Clonality: | Monoclonal Antibody |
| SwissProt ID: | O00151 |
| Immunogen: | A synthesized peptide derived from human PDLIM1 |
| Gene ID: | 9124 |

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: PDLIM1
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 36 kDa, Observed MW: 36 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

| Antibody Dilution Ratio: | Application | Antibody Dilution Ratio |
|---------------------------------|--------------------|--------------------------------|
| | WB | 1:1000-1:2000 |
| | IF | 1:50-1:200 |
| | FC | 1:50-1:100 |
| | IP | 1:50 |

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.