

Anti-Tat SF1 [R09-4L3] Monoclonal Antibody

AGMB03285

Description

This Anti-Tat SF1 [R09-4L3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03285
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	HTATSF1, HTSF1, TAT SF1
Applications:	WB IHC-P ICC/IF IP
Research Area:	Microbiology
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse
Clone:	R09-4L3
Clonality:	Monoclonal Antibody
SwissProt ID:	O43719
Immunogen:	A synthesized peptide derived from human HTSF1
Gene ID:	27336
Gene Name:	HTATSF1

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 86 kDa, Observed MW: 140 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC	1:50-1:100
	IF	1:50-1:200
	IP	1:50

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.