

Anti-IL-32 [R09-1M3] Monoclonal Antibody

AGMB03293

Description

This Anti-IL-32 [R09-1M3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

| | |
|-----------------------|---|
| SKU: | AGMB03293 |
| Contents: | 20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml) |
| Synonyms: | IL-32, IL 32alpha, IL 32beta, IL 32delta, IL 32gamma, IL32, NK4, TAIF, TAIFa, TAIFb, TAIFc, TAIFd |
| Applications: | WB |
| Research Area: | Immunology |
| Form: | Liquid |

Antibody Data

| | |
|----------------------|---|
| Reactivity: | Human |
| Clone: | R09-1M3 |
| Clonality: | Monoclonal Antibody |
| SwissProt ID: | P24001 |
| Immunogen: | A synthesized peptide derived from human IL32 |
| Gene ID: | 9235 |

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: IL32
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 27 kDa, Observed MW: 23 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

| Antibody Dilution Ratio: | Application | Antibody Dilution Ratio |
|---------------------------------|--------------------|--------------------------------|
| | WB | 1:500-1:1000 |

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.