

Anti-TAF1 [R09-7P2] Monoclonal Antibody

AGMB03302

Description

This Anti-TAF1 [R09-7P2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03302
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	BA2R, CCG1, CCGS, Cell cycle gene 1 protein, DYT3, KAT4, N TAF1, NSCL2, OF, P250, TAF(II)250, TAF1, TAF1-Specific, TAF2A, TAFII 250, TAFII250, TBP associated factor 250 kDa
Applications:	WB
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R09-7P2
Clonality:	Monoclonal Antibody
SwissProt ID:	P21675
Immunogen:	A synthesized peptide derived from human TAF1
Gene ID:	6872

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: TAF1
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 213 kDa, Observed MW: 300 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

	Application	Antibody Dilution Ratio
Antibody Dilution Ratio:	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.