

Anti-CBR1 [R06-6L5] Monoclonal Antibody

AGMB03313

Description

This Anti-CBR1 [R06-6L5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03313
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	15 hydroxyprostaglandin dehydrogenase [NADP], 15-hydroxyprostaglandin dehydrogenase [NADP], Carbonyl reductase [NADPH] 1, CBR 1, CBR1, CBR1_HUMAN, CRN, NADPH dependent carbonyl reductase 1, NADPH-dependent carbonyl reductase 1, Prostaglandin 9 ketoreductase, Prostaglandin 9-ketoreductase, Prostaglandin E(2) 9 reductase, Prostaglandin-E(2) 9-reductase, SDR21C1.
Applications:	
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R06-6L5
Clonality:	Monoclonal Antibody
SwissProt ID:	P16152

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Immunogen:	A synthesized peptide derived from human CBR1
Gene ID:	873
Gene Name:	CBR1
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 30 kDa, Observed MW: 30 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IP	1:50

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.