

Anti-ARID1A [R07-8X2] Monoclonal Antibody

AGMB03322

Description

This Anti-ARID1A [R07-8X2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03322
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Actin-dependent regulator of chromatin subfamily F member 1, ARI1A_HUMAN, ARID domain containing protein 1A,
Applications:	IHC-P ICC/IF
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R07-8X2
Clonality:	Monoclonal Antibody
SwissProt ID:	O14497
Immunogen:	A synthesized peptide derived from human ARID1A
Gene ID:	8289

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: ARID1A

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 242 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:100-1:200
	IF	1:50-1:200
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	