

# Anti-Lambda Light chain [R51-1S-3] Monoclonal Antibody

AGMB03346

## Description

---

This Anti-Lambda Light chain [R51-1S-3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

<b>SKU:</b>	AGMB03346
<b>Contents:</b>	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	Bence Jones protein, BJP, Constant region of lambda light chains, Ig lambda chain C regions, IGLC, IGLC1, IGLC2, IGLC3, IGLC6, IGLV, Immunoglobulin lambda constant 1, Mcg marker, Paraprotein,
<b>Applications:</b>	<b>WB</b> <b>IHC-P</b> <b>ICC/IF</b> <b>FC</b>
<b>Research Area:</b>	Immunology
<b>Form:</b>	Liquid

## Antibody Data

---

<b>Reactivity:</b>	Human
<b>Clone:</b>	R51-1S-3
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P0CG04
<b>Immunogen:</b>	A synthesized peptide derived from human Lambda Light chain
<b>Gene ID:</b>	28820

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** IGLV1-51  
**Host Species:** Rabbit  
**Isotype:** IgG  
**Purification:** Affinity Chromatography  
**Conjugated:** Unconjugated  
**Modification:** Unmodified  
**Molecular Weight:** Calculated MW: 12 kDa, Observed MW: 25kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

<b>Antibody Dilution Ratio:</b>	<b>Application</b>	<b>Antibody Dilution Ratio</b>
	WB	1:1000-1:2000
	IHC	1:100-1:200
	IF	1:50-1:200
	FC	1:50

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.