

# Anti-Sodium/Hydrogen Exchanger 1/NHE-1 [R86-4H-8] Monoclonal Antibody

AGMB03354

## Description

---

This Anti-Sodium/Hydrogen Exchanger 1/NHE-1 [R86-4H-8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

<b>SKU:</b>	AGMB03354
<b>Contents:</b>	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	SLC9A1, APNH, APNH1, NHE1, LIKNS, NHE-1, PPP1R143, Sodium/hydrogen exchanger 1
<b>Applications:</b>	<b>WB</b>
<b>Research Area:</b>	Signal Transduction
<b>Form:</b>	Liquid

## Antibody Data

---

<b>Reactivity:</b>	Human
<b>Clone:</b>	R86-4H-8
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P19634
<b>Immunogen:</b>	A synthesized peptide derived from human Sodium / Hydrogen Exchanger 1

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene ID:** 6548  
**Gene Name:** SLC9A1  
**Host Species:** Rabbit  
**Isotype:** IgG  
**Purification:** Affinity Chromatography  
**Conjugated:** Unconjugated  
**Modification:** Unmodified  
**Molecular Weight:** Calculated MW: 91 kDa, Observed MW: 130kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.