

Anti-SAHH [R35-5L-9] Monoclonal Antibody

AGMB03363

Description

This Anti-SAHH [R35-5L-9] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03363			
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)			
Synonyms:	Adenosylhomocysteinase,	AdoHcyase,	ahcY,	SAHH,
	Adenosylhomocysteinase,	AdoHcyase,	ahcY,	
Applications:				
Research Area:	Signal Transduction			
Form:	Liquid			

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R35-5L-9
Clonality:	Monoclonal Antibody
SwissProt ID:	P23526
Immunogen:	A synthesized peptide derived from AHCY
Gene ID:	191

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: AHCY

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 48 kDa, Observed MW: 47 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.