

Anti-Erlin-2 [R42-7O-2] Monoclonal Antibody

AGMB03384

Description

This Anti-Erlin-2 [R42-7O-2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03384
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	ERLIN2, NET32, SPFH2, C8orf2, Endoplasmic reticulum lipid raft-associated protein 2, Stomatin-prohibitin-flotillin-HflC/K domain-containing protein 2, SPFH domain-containing protein 2
Applications:	WB IHC-P ICC/IF
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse
Clone:	R42-7O-2
Clonality:	Monoclonal Antibody
SwissProt ID:	O94905
Immunogen:	A synthesized peptide derived from human Erlin-2
Gene ID:	11160

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: ERLIN2
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Chromatography
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 38 kDa, Observed MW: 38kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC	1:100-1:200
	IF	1:50-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.