

Anti-ADH4 [R94-5J-6] Monoclonal Antibody

AGMB03430

Description

This Anti-ADH4 [R94-5J-6] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03430
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Alcohol dehydrogenase 4, Alcohol dehydrogenase class II pi chain, Ac1002, ADH 4, ADH-1,
Applications:	WB IHC-P ICC/IF
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R94-5J-6
Clonality:	Monoclonal Antibody
SwissProt ID:	P08319
Immunogen:	A synthesized peptide derived from human ADH4
Gene ID:	127

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: ADH4
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Chromatography
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW:40kDa, Observed MW:37kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC	1:100-1:200
	IF	1:50-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.