

Anti-UQCRH [R42-3F-2] Monoclonal Antibody

AGMB03431

Description

This Anti-UQCRH [R42-3F-2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03431
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Complex III subunit 6, Complex III subunit VIII, Cytochrome b c1 complex subunit 6, QCR6, UQCR8, MC3DN11
Applications:	WB IHC-P
Research Area:	Tags & Cell Markers
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R42-3F-2
Clonality:	Monoclonal Antibody
SwissProt ID:	P07919
Immunogen:	A synthesized peptide derived from human UQCRH
Gene ID:	7388

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	UQCRH
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Chromatography
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW:11kDa,Observed MW:11kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC	1:100-1:200
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	