

# Anti-Somatostatin Receptor 1/SSTR1 [R95-4Q-5] Monoclonal Antibody AGMB03432

## Description

---

This Anti-Somatostatin Receptor 1/SSTR1 [R95-4Q-5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

<b>SKU:</b>	AGMB03432
<b>Contents:</b>	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	SS1R, SST1, SS1-R, SRIF-2, SS-1-R, SSTR1
<b>Applications:</b>	<b>WB</b> <b>ICC/IF</b>
<b>Research Area:</b>	Neuroscience
<b>Form:</b>	Liquid

## Antibody Data

---

<b>Reactivity:</b>	Human
<b>Clone:</b>	R95-4Q-5
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P30872
<b>Immunogen:</b>	A synthesized peptide derived from human Somatostatin Receptor 1
<b>Gene ID:</b>	6751
<b>Gene Name:</b>	SSTR1

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Host Species:** Rabbit

**Isotype:** IgG

**Purification:** Affinity Chromatography

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** Calculated MW:43kDa, Observed MW:60kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IF	1:50-1:200

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.