

Anti-TOE1 [R11-5F-2] Monoclonal Antibody

AGMB03437

Description

This Anti-TOE1 [R11-5F-2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

| | |
|-----------------------|---|
| SKU: | AGMB03437 |
| Contents: | 20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml) |
| Synonyms: | Toe1, Target of EGR1 protein 1 |
| Applications: | WB |
| Research Area: | Cell Biology |
| Form: | Liquid |

Antibody Data

| | |
|----------------------|---|
| Reactivity: | Human |
| Clone: | R11-5F-2 |
| Clonality: | Monoclonal Antibody |
| SwissProt ID: | Q96GM8 |
| Immunogen: | A synthesized peptide derived from human TOE1 |
| Gene ID: | 114034 |
| Gene Name: | TOE1 |

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 57 kDa, Observed MW: 60kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

| Antibody Dilution Ratio: | Application | Antibody Dilution Ratio |
|--------------------------|-------------|-------------------------|
| | WB | 1:500-1:1000 |

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.