

Anti-CD32A+CD32B+CD32C [R29-3T-5] Monoclonal Antibody

AGMB03448

Description

This Anti-CD32A+CD32B+CD32C [R29-3T-5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03448
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CD32, CD32A, CD32B, CD32C, CDw32, FCG2, FcGR, FCGR2, FCGR2A, FCGR2A1, FCGR2B, FCGR2C, Fcgr3, Fcgr3a, FcgrII, FcRII a, FCRII, FcRII b, FcRII c, FcRIIC, IGFR2,
Applications:	WB IHC-P
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R29-3T-5
Clonality:	Monoclonal Antibody
SwissProt ID:	P12318/P31994/P31995
Immunogen:	A synthesized peptide derived from human CD32A + CD32B + CD32C
Gene ID:	2212/2213/9103

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: FCGR2A/FCGR2B/FCGR2C

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW:34/35/36kDa, Observed MW:37-43kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC	1:100-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.