

Anti-SOCS7 [R27-4V-7] Monoclonal Antibody

AGMB03455

Description

This Anti-SOCS7 [R27-4V-7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03455
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Ash and phospholipase C gamma-binding protein, NAP4, Nck, Nck Ash and phospholipase C binding protein, Nck associated protein 4, NCKAP4, SH2 domain containing SOCS box protein, SOCS4, SOCS6, Socs7, Suppressor of cytokine signaling 7,
Applications:	WB IHC-P
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R27-4V-7
Clonality:	Monoclonal Antibody
SwissProt ID:	O14512
Immunogen:	A synthesized peptide derived from human SOCS7

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID: 30837
Gene Name: SOCS7
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Chromatography
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW:63kDa, Observed MW:63kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC	1:100-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.