

Anti-RIP2 [R87-7G-8] Monoclonal Antibody

AGMB03458

Description

This Anti-RIP2 [R87-7G-8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03458
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CARD carrying kinase, CARD3, CARDIAK, CCK, CLARP kinase, GIG30, Growth inhibiting gene 30, Receptor interacting protein 2, RICK, RIP 2, Ripk2, TNFRSF
Applications:	WB
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R87-7G-8
Clonality:	Monoclonal Antibody
SwissProt ID:	O43353
Immunogen:	A synthesized peptide derived from human RIP2
Gene ID:	8767

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: RIPK2
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Chromatography
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW:61kDa, Observed MW:61kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

	Application	Antibody Dilution Ratio
Antibody Dilution Ratio:	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.