

Anti-DcR1/CD263 [R67-7H-8] Monoclonal Antibody

AGMB03471

Description

This Anti-DcR1/CD263 [R67-7H-8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03471
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Antagonist decoy receptor for TRAIL/Apo 2L, Antagonist decoy receptor for TRAIL/Apo-2L, CD263,
Applications:	WB ICC/IF
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R67-7H-8
Clonality:	Monoclonal Antibody
SwissProt ID:	O14798
Immunogen:	A synthesized peptide derived from human DcR1 / CD263
Gene ID:	8794

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: TNFRSF10C

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW:27kDa, Observed MW:27kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IF	1:50-1:200
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	