

Anti-Annexin-7/ANXA7 [R03-3D-3] Monoclonal Antibody

AGMB03517

Description

This Anti-Annexin-7/ANXA7 [R03-3D-3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03517
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Annexin VII, SNX, ANX7, SYNEXIN, ANXA7, Annexin A7, Annexin VII, Annexin-7,
Applications:	IHC-P
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R03-3D-3
Clonality:	Monoclonal Antibody
SwissProt ID:	P20073
Immunogen:	Recombinant protein of human Annexin VII
Gene ID:	310

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: ANXA7
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:100-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.