

# Anti-ErbB3 / HER3 [R06-3O-3] Monoclonal Antibody

AGMB03602

## Description

---

This Anti-ErbB3 / HER3 [R06-3O-3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

|                       |  |
|-----------------------|--|
| <b>SKU:</b>           | AGMB03602  |
| <b>Contents:</b>      | 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml)  |
| <b>Synonyms:</b>      | ERBB3, HER3, Receptor tyrosine-protein kinase erbB-3, Proto-oncogene-like protein c-ErbB-3, Tyrosine kinase-type cell surface receptor HER3, c erbB 3, c erbB3, Erb b2 receptor tyrosine kinase 3, |
| <b>Applications:</b>  | <b>IHC-P</b>   |
| <b>Research Area:</b> | Signal Transduction  |
| <b>Form:</b>          | Liquid   |

## Antibody Data

---

|                      |  |
|----------------------|--|
| <b>Reactivity:</b>   | Human  |
| <b>Clone:</b>        | R06-3O-3                                       |
| <b>Clonality:</b>    | Monoclonal Antibody                            |
| <b>SwissProt ID:</b> | P21860   |
| <b>Immunogen:</b>    | A synthesized peptide derived from human ERBB3 |
| <b>Gene ID:</b>      | 2065   |

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** ERBB3  
**Host Species:** Rabbit  
**Isotype:** IgG  
**Purification:** Affinity Purified  
**Conjugated:** Unconjugated  
**Modification:** Unmodified  
**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

| <b>Antibody Dilution Ratio:</b> | <b>Application</b> | <b>Antibody Dilution Ratio</b> |
|---------------------------------|--------------------|--------------------------------|
|                                 | IHC                | 1:100-1:500                    |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.