

Anti-FOXM1 [R09-2P-3] Monoclonal Antibody

AGMB03614

Description

This Anti-FOXM1 [R09-2P-3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03614
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	MPP2, HFH11, HNF-3, INS-1, MPP-2, PIG29, FKHL16, FOXM1A, FOXM1B, FOXM1C, HFH-11, TRIDENT, MPHOSPH2
Applications:	IHC-P
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R09-2P-3
Clonality:	Monoclonal Antibody
SwissProt ID:	Q08050
Immunogen:	Recombinant protein of human FOXM1
Gene ID:	2305

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: FOXM1
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:100-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.