

Anti-NKX6.1 [R09-9I-2] Monoclonal Antibody

AGMB03697

Description

This Anti-NKX6.1 [R09-9I-2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03697
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	NKX6-1, NKX6A, Homeobox protein Nkx-6.1, Homeobox protein NK-6 homolog A, Homeobox protein NK-6 homolog A, Homeobox protein Nkx-6.1, NK homeo box, family 6, member A,
Applications:	IHC-P
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R09-9I-2
Clonality:	Monoclonal Antibody
SwissProt ID:	P78426
Immunogen:	Recombinant protein of human NKX6.1
Gene ID:	4825

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: NKX6-1
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

	Application	Antibody Dilution Ratio
Antibody Dilution Ratio:	IHC-P	1:1000-1:2000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.