

Anti-RBX1 [R07-2U-1] Monoclonal Antibody

AGMB03726

Description

This Anti-RBX1 [R07-2U-1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03726
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	ROC1, RNF75, BA554C12.1, RBX1, BA554C12.1, E3 ubiquitin-protein ligase RBX1, FLJ60363,
Applications:	IHC-P
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R07-2U-1
Clonality:	Monoclonal Antibody
SwissProt ID:	P62877
Immunogen:	A synthesized peptide derived from human ROC1
Gene ID:	9978

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: RBX1
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC-P	1:2000-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.