

# Anti-SLC7A5 [R05-7J-7] Monoclonal Antibody

AGMB03744

## Description

---

This Anti-SLC7A5 [R05-7J-7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

<b>SKU:</b>	AGMB03744
<b>Contents:</b>	50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	SLC7A5, CD98LC, LAT1, MPE16, Large neutral amino acids transporter small subunit 1, 4F2 light chain (4F2 LC, 4F2LC), CD98 light chain, Integral membrane protein E16, Solute carrier family 7 member 5
<b>Applications:</b>	<b>IHC-P</b>
<b>Research Area:</b>	Signal Transduction
<b>Form:</b>	Liquid

## Antibody Data

---

<b>Reactivity:</b>	Human
<b>Clone:</b>	R05-7J-7
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	Q01650
<b>Immunogen:</b>	A synthesized peptide derived from human SLC7A5
<b>Gene ID:</b>	8140

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** SLC7A5  
**Host Species:** Rabbit  
**Isotype:** IgG  
**Purification:** Affinity Purified  
**Conjugated:** Unconjugated  
**Modification:** Unmodified  
**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

	Application	Antibody Dilution Ratio
<b>Antibody Dilution Ratio:</b>	IHC	1:100-1:200

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.