

Anti-Somatostatin Receptor 2 [R01-1P-5] Monoclonal Antibody

AGMB03746

Description

This Anti-Somatostatin Receptor 2 [R01-1P-5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03746
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	SSTR2, somatostatin receptor 2, SS 2 R, SS2R, SRIF 1, SS2 R
Applications:	IHC-P
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R01-1P-5
Clonality:	Monoclonal Antibody
SwissProt ID:	P30874
Immunogen:	A synthesized peptide derived from human Somatostatin Receptor 2
Gene ID:	6752
Gene Name:	SSTR2

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC-P	1:500-1:2000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.