

Anti-TRP1 [R01-1N-9] Monoclonal Antibody

AGMB03758

Description

This Anti-TRP1 [R01-1N-9] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03758
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CAS2, CATB, GP75, OCA3, TRP1, TYRP, TYRP1, TYRRP, 5, 5,6 dihydroxyindole 2 carboxylic acid oxidase, 6-dihydroxyindole-2-carboxylic acid oxidase,
Applications:	IHC-P
Research Area:	Tags & Cell Markers
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R01-1N-9
Clonality:	Monoclonal Antibody
SwissProt ID:	P17643
Immunogen:	Recombinant protein of human TRP1
Gene ID:	7306

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: TYRP1
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC-P	1:2000-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.