

## Anti-Bcl 6 [R06-5S-4] Monoclonal Antibody

AGMB03772

### Description

---

This Anti-Bcl 6 [R06-5S-4] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |                                                                                                                                                     |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>SKU:</b>           | AGMB03772                                                                                                                                           |
| <b>Contents:</b>      | 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml)                                                                                                       |
| <b>Synonyms:</b>      | BCL6, BCL5, LAZ3, ZBTB27, ZNF51, B-cell lymphoma 6 protein, Protein LAZ-3, Zinc finger and BTB domain-containing protein 27, Zinc finger protein 51 |
| <b>Applications:</b>  | <b>IHC-P</b>                                                                                                                                        |
| <b>Research Area:</b> | Epigenetics                                                                                                                                         |
| <b>Form:</b>          | Liquid                                                                                                                                              |

### Antibody Data

---

|                      |                                    |
|----------------------|------------------------------------|
| <b>Reactivity:</b>   | Human                              |
| <b>Clone:</b>        | R06-5S-4                           |
| <b>Clonality:</b>    | Monoclonal Antibody                |
| <b>SwissProt ID:</b> | P41182                             |
| <b>Immunogen:</b>    | Recombinant protein of human BCL-6 |
| <b>Gene ID:</b>      | 604                                |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** BCL6

**Host Species:** Rabbit

**Isotype:** IgG

**Purification:** Affinity Purified

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

| <b>Antibody Dilution Ratio:</b> | <b>Application</b> | <b>Antibody Dilution Ratio</b> |
|---------------------------------|--------------------|--------------------------------|
|                                 | IHC-P              | 1:200-1:2000                   |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.