

Anti-CDK1 [R03-7B-5] Monoclonal Antibody

AGMB03778

Description

This Anti-CDK1 [R03-7B-5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03778
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CDK1, CDC2, CDC28A, CDKN1, P34CDC2, Cyclin-dependent kinase 1, CDK1, Cell division control protein 2 homolog, Cell division protein kinase 1, p34 protein kinase
Applications:	ICC/IF FC
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R03-7B-5
Clonality:	Monoclonal Antibody
SwissProt ID:	P06493
Immunogen:	Peptide
Gene ID:	983

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: CDK1
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at -20°C. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IF	1:50-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.