

Anti-E Cadherin [R08-3J-8] Monoclonal Antibody

AGMB03786

Description

This Anti-E Cadherin [R08-3J-8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03786
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CDH1, CDHE, UVO, Cadherin-1, CAM 120/80, Epithelial cadherin, E-cadherin, Uvomorulin, CD antigen CD324
Applications:	WB
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R08-3J-8
Clonality:	Monoclonal Antibody
SwissProt ID:	P12830
Immunogen:	Peptide
Gene ID:	999

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: CDH1
Host Species: Rabbit
Isotype: IgG
Purification: Protein A
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 97 kDa, Observed MW:90-120 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.