

Anti-STAT1 alpha [R09-8C-1] Monoclonal Antibody

AGMB03791

Description

This Anti-STAT1 alpha [R09-8C-1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03791
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Stat1, Signal transducer and activator of transcription 1
Applications:	WB IHC-P IP
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R09-8C-1
Clonality:	Monoclonal Antibody
SwissProt ID:	P42225
Immunogen:	Peptide
Gene ID:	20846
Gene Name:	Stat1

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species:	Rabbit
Isotype:	IgG
Purification:	Protein A
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 87 kDa, Observed MW:91 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:2000
	IHC-P	1:50-1:100
	IP	1:25-1:50
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	