

Anti-CD42b [R07-5I-7] Monoclonal Antibody

AGMB03798

Description

This Anti-CD42b [R07-5I-7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03798
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Antigen CD42b alpha, BDPLT1, BDPLT3, BSS, CD42b alpha, DBPLT3, Glycoprotein Ibalpha, GP Ib alpha, GP1B, GP1BA, GPIb alpha, GPIbA, VWDP
Applications:	IHC-P
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R07-5I-7
Clonality:	Monoclonal Antibody
SwissProt ID:	P07359
Immunogen:	Peptide
Gene ID:	2811

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: GP1BA
Host Species: Rabbit
Isotype: IgG
Purification: Protein A
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 72 kDa, Observed MW:120 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC-P	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.