

Anti-CD58 [R07-5X-4] Monoclonal Antibody

AGMB03800

Description

This Anti-CD58 [R07-5X-4] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03800
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	AG3, CD58, Surface glycoprotein LFA-3, CD58 antigen (Lymphocyte function-associated antigen 3),
Applications:	WB ICC/IF FC
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R07-5X-4
Clonality:	Monoclonal Antibody
SwissProt ID:	P19256
Immunogen:	Peptide
Gene ID:	965

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: CD58
Host Species: Rabbit
Isotype: IgG
Purification: Protein A
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 28 kDa, Observed MW:65-70 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	ICC/IF	1:50
	FC	1:50-1:100
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	