

Anti-CD23 [R07-3V-8] Monoclonal Antibody

AGMB03809

Description

This Anti-CD23 [R07-3V-8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03809
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	FCER2, CD23A, CLEC4J, FCE2, IGEBF, Low affinity immunoglobulin epsilon Fc receptor, BLAST-2, C-type lectin domain family 4 member J, Fc-epsilon-RII, Immunoglobulin E-binding factor, Lymphocyte IgE receptor, CD23
Applications:	IHC-P ICC/IF FC
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R07-3V-8
Clonality:	Monoclonal Antibody
SwissProt ID:	P06734
Immunogen:	Peptide
Gene ID:	2208

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: FCER2

Host Species: Rabbit

Isotype: IgG

Purification: Protein A

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC-P	1:50-1:100
	ICC/IF	1:50-1:100
	FC	1:50-1:100
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	