

Anti-ALOX15 [R06-4V-7] Monoclonal Antibody

AGMB03814

Description

This Anti-ALOX15 [R06-4V-7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03814
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	ALOX15, LOG15, LOG15, 12-LOX, 15-LOX, 15-LOX-1, Polyunsaturated fatty acid lipoxygenase ALOX15
Applications:	IHC-P ICC/IF FC
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R06-4V-7
Clonality:	Monoclonal Antibody
SwissProt ID:	P16050
Immunogen:	Peptide
Gene ID:	246

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: ALOX15

Host Species: Rabbit

Isotype: IgG

Purification: Protein A

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC-P	1:50-1:100
	ICC/IF	1:50
	FC	1:50-1:100
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	