

Anti-SFTP A1 + SFTP A2 [R09-7P-8] Monoclonal Antibody

AGMB03828

Description

This Anti-SFTP A1 + SFTP A2 [R09-7P-8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03828
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Collectin-4, COLEC4, Collectin-5, COLEC5, PSAP, SFTP1, SFTP A, SFTP A1B, SFTP A2B
Applications:	
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R09-7P-8
Clonality:	Monoclonal Antibody
SwissProt ID:	Q8IWL2, Q8IWL1
Immunogen:	Peptide
Gene ID:	653509;729238

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: SFTPA1, SFTPA2

Host Species: Rabbit

Isotype: IgG

Purification: Protein A

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 27 kDa, Observed MW:27 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.