

Anti-Sialoadhesin/CD169 [R07-2E-1] Monoclonal Antibody

AGMB03831

Description

This Anti-Sialoadhesin/CD169 [R07-2E-1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03831
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	SIGLEC1, SN, SIGLEC-1, dJ1009E24.1, CD169, CD169 antigen, dJ1009E24.1
Applications:	IHC-P ICC/IF FC
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R07-2E-1
Clonality:	Monoclonal Antibody
SwissProt ID:	Q9BZZ2
Immunogen:	Peptide
Gene ID:	6614
Gene Name:	SIGLEC1

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Rabbit
Isotype: IgG
Purification: Protein A
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC-P	1:50-1:100
	ICC/IF	1:50
	FC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.