

Anti-Cytokeratin 20 [R08-8V-8] Monoclonal Antibody

AGMB03839

Description

This Anti-Cytokeratin 20 [R08-8V-8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03839
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CD20, CK-20, CK20, cytokeratin 20, Cytokeratin-20, K1C20, K20, keratin 20, Keratin, type I cytoskeletal 20, Keratin-20, KRT20, KRT21, MGC35423, Protein IT
Applications:	WB IHC-P ICC/IF FC
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R08-8V-8
Clonality:	Monoclonal Antibody
SwissProt ID:	P35900
Immunogen:	Peptide
Gene ID:	54474

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: KRT20

Host Species: Rabbit

Isotype: IgG

Purification: Protein A

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 48 kDa, Observed MW:48 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC-P	1:50-1:100
	ICC/IF	1:50
	FC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.