

Anti-BAP31 [R06-7O-4] Monoclonal Antibody

AGMB03840

Description

This Anti-BAP31 [R06-7O-4] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03840
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	BCAP31, BAP31, DXS1357E, B-cell receptor-associated protein 31, BCR-associated protein 31, Bap31, 6C6-AG tumor-associated antigen, Protein CDM, p28
Applications:	WB
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R06-7O-4
Clonality:	Monoclonal Antibody
SwissProt ID:	P51572
Immunogen:	Peptide
Gene ID:	10134

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: BCAP31

Host Species: Rabbit

Isotype: IgG

Purification: Protein G

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 28 kDa, Observed MW:28 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.