

# Anti-PPAR gamma [R06-7T-7] Monoclonal Antibody

AGMB03859

## Description

---

This Anti-PPAR gamma [R06-7T-7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

**SKU:** AGMB03859

**Contents:** 20µl, 50µl, 100µl  
Bradford Reagent: 1 vial (2ml)

**Synonyms:** CIMT1, GLM1, NR1C3, Nuclear receptor subfamily 1 group C member 3, Peroxisome proliferator activated nuclear receptor gamma variant 1, Peroxisome proliferator activated receptor gamma 1, Peroxisome Proliferator Activated Receptor gamma, Peroxisome proliferator-activated receptor gamma, PPAR gamma, PPAR-gamma, PPARG, PPARG\_HUMAN, PPARG1, PPARG2, PPARGgamma.

**Applications:** **WB** **ICC/IF** **FC**

**Research Area:** Epigenetics

**Form:** Liquid

## Antibody Data

---

**Reactivity:** Human

**Clone:** R06-7T-7

**Clonality:** Monoclonal Antibody

**SwissProt ID:** P37231

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Immunogen:** Peptide

**Gene ID:** 5468

**Gene Name:** PPARG

**Host Species:** Rabbit

**Isotype:** IgG

**Purification:** Affinity Purified

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** Calculated MW:58 kDa,Observed MW: 58 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

<b>Antibody Dilution Ratio:</b>	<b>Application</b>	<b>Antibody Dilution Ratio</b>
	WB	1:500-1:1000
	ICC/IF	1:100-1:200
	FC	1:200-1:1000

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.