

Anti-MKK7 [R03-2N-8] Monoclonal Antibody

AGMB03868

Description

This Anti-MKK7 [R03-2N-8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU: AGMB03868

Contents: 20µl, 50µl, 100µl
Bradford Reagent: 1 vial (2ml)

Synonyms: MEK, MKK7, JNKK2, MEK 7, MAPKK7, PRKMK7, SAPKK4, SAPKK-4, K7, c-Jun N-terminal kinase kinase 2, Dual specificity mitogen activated protein kinase kinase 7, Dual specificity mitogen-activated protein kinase kinase 7,

Applications: 

Research Area: Immunology

Form: Liquid

Antibody Data

Reactivity: Human, Mouse, Rat

Clone: R03-2N-8

Clonality: Monoclonal Antibody

SwissProt ID: O14733

Immunogen: Peptide

Gene ID: 5609

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: MAP2K7
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW:48 kDa, Observed MW: 48 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	IHC-P	1:50-1:100
	FC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.